

porodo

| CHOOSE THE BEST

Accupressure Hand Massager
With Six Mode

SKU:PD-APHNDMS-GY



Before installing and using the product, please carefully read this User Manual to guarantee correct usage.

(Please keep this instruction manual securely for future use.)

Charge the device upon purchase using a 5V, 1A (or 2A) cell phone charger. The power indicator will flash during charging and will turn off once charging is complete, indicating full battery.


Suitable for

- Individuals who use computers for extended periods.
- Those with poor palm blood circulation.
- Women who wish to maintain hand beauty.

Product Features

- *Ergonomic design tailored to hand's acupressure points.
- *Pneumatic segmentation for massage.
- *Comprehensive palm acupuncture massage.
- *Heat function for comfort in cold conditions (optional).
- *Microseismic massage for enhanced experience (optional).

Control Panel



Features



Air Pressure

Auto

Auto



Heating



Vibration


Air pressure massage feature.

Automated massage programs.

Heat therapy option.

Vibration massage feature.

Components



How to Use the Control Panel



Power Button

Press for 2 seconds to turn on/off.

The device runs for 15 minutes by default.



Intensity Button

Has 6 intensity levels, defaulting to level 3 (L3).

(L1, L2, L3, L4, L5, L6) (L1: Lowest, L2: Highest)

How to Use the Control Panel



Mode Key

Provides 6 massage modes, starting at mode P1. (P1, P2, P3, P4, P5, P6)



Vibration Button

Toggle vibration feature with 2 modes available.



Heating Key

Activates heat function, noticeable increase in temperature after 5 minutes.

Troubleshooting

1. If the device won't start, ensure it's charged or press the power button for 2 seconds.
2. Normal sounds include air pump buzzing and solenoid valve ticking, indicating no fault.

Specifications

Product Name	Hand Massager With Six Mode	Power	5W
Model	PD-APHNDMS	Operational Time	15 minutes
Voltage	DC 3.7V	Weight	1.1 kg